



VACC 201 Tinkering and Mentoring

Name of the Product

TriSense-Heal Patch

A multi-biomarker responsive smart wound dressing for early infection detection



Early Biochemical Warning Before Visible Infection

Only plant/bio-sourced biodegradable materials used



Mentor : Dr. Yesha Master

Name	Enrollment Number	Course
Aditi Singh	25SC04001	BSc Microbiology
Ali Asger Sakir Husain Anandwala	25SC04002	BSc Microbiology
Anup Rajendra Tanwar	25SC04003	BSc Microbiology
Bapat Vrinda Kedar	25SC04004	BSc Microbiology
Chaudhari Dhruvikumari Shaileshbhai	25SC04005	BSc Microbiology
Dixit Devanshi Hemkant	25SC04006	BSc Microbiology

Objectives

Primary Objectives

To develop a **biodegradable bandage system** using agar-based hydrogel and natural materials

To integrate **multi-parameter sensing**, including:

- pH changes (infection indicator)
- Protease activity (tissue degradation marker)
- Oxidative stress (healing/inflammation marker)
- Moisture level (wound environment)
- Temperature variation (inflammation/infection indicator)

Secondary Objectives

To utilize **natural bioactive compounds**:

- Anthocyanin (pH sensing)
- Bromelain (proteolytic activity & wound cleaning)
- Curcumin (antimicrobial & antioxidant)

To design a **hybrid architecture**:

- Functional hydrogel layer (healing)
- Surface sensor zones (visual detection)

To ensure:

- Cost-effectiveness**
- Ease of use (visual color-based detection)**
- Scalability for mass production**

Applied Objective


To create a **prototype capable of early infection detection**, reducing dependence on delayed visual diagnosis and enabling timely medical intervention.

Problem Statement

Conventional bandages are passive, non-biodegradable, and lack real-time monitoring, leading to delayed infection detection. This creates a need for a biodegradable smart bandage capable of monitoring multiple wound parameters simultaneously.

CONCEPT OVERVIEW (GRAPHICAL REPRESENTATION)

The TriSense-Heal Patch



ARCHITECTURE

Low-Cost, Natural Biomaterial Sourcing

Replacing expensive synthetic lab chemicals with highly effective, biodegradable natural equivalents.

Bioactive Components



Turmeric (Curcumin)
- Antioxidant properties and ROS redox response.



Pineapple Extract (Bromelain)
- Natural enzyme integration for protease interactions.

Structural Base



Gelatin / Alginate (100–200 g)
- Hydrogel base for sensing zones (Approx. Cost: 800 INR).



Sterile Cellulose Gauze
- Soft, absorbent contact layer.

Sensor Indicators



Red Cabbage Extract (Anthocyanin)
- Natural pH-sensitive colorimetric indicator.



Resazurin / Bromothymol Blue
- Supplemental oxidative and pH validation dyes.

Laboratory Prototype Development

Translating accessible raw biomaterials into a precisely engineered, functional diagnostic dressing.



1. Extraction: Beakers containing vibrant extracted natural anthocyanin dyes.



2. Base Formation: Liquid hydrogel poured into precise acrylic molds.



3. Sensor Integration: Micro-dispensing the circular sensor zones.




4. Final Assembly: The fully assembled prototype patch in its dry state.

PROTOTYPE DEVELOPMENT



PROTOTYPE




PROTEASE LAYER



FOIL PACKAGING

ILLUSTRATION ON PACKAGING

TriSense-Heal™ Patch



Response interpretation


Green = Normal Yellow = Warning Red = Severe

OBSERVATION TABLE OF VISIBLE RESPONSE

Sensor		Trigger	Before	After
1	pH Variation	Wound shifts to alkaline	 Purple	 Green/Yellow
2	Protease Activity	Bacterial enzymes degrade gelatin	 Contained Dye	 Matrix Breakdown & Spread
3	Oxidative Stress	Elevated ROS from inflammation	 Baseline	 Intensity Change via Redox
4	Moisture	Excessive exudate buildup	 Dry Strip	 Physical Swelling
5	Temperature	Localized inflammatory heat	 Neutral	 Thermochromic Shift

CHALLENGES AND IMPROVEMENTS

Hydrated Gel Shelf-Life Limitations



Dry-State Activation

The patch remains inert during storage and activates only upon contact with wound exudate.

Natural Biomolecule Degradation



Micro-Encapsulation

Shielding active ingredients like curcumin to prevent premature fading and ensure long-term molecular stability.

Environmental Contamination Risk



Advanced Packaging

Utilizing medical-grade vacuum sealing to ensure complete sterility and baseline stability until application.




GSFC
UNIVERSITY
EDUCATION RE-ENVISIONED




GUITAR
COUNCIL

WOUND CARE – SMART WOUND CARE




SCAN TO TRACK YOUR WOUND AND
PREDICT INFECTION RATE



The Future of Diagnostic Wound Care

A scalable, single-use consumable that shifts wound care from passive observation to proactive medical intervention.

The Value Triangle



Development Roadmap



...THANK YOU